



FACE MASK DISPOSABLE

~~\$2.70~~ **\$1.79**

Certified Disposable Face Mask
Minimum Order: 20 Masks

FREE DELIVERY

SKU: DCC-95

Categories: [Bioide](#), [FaceProtect](#), [Masks](#),
[PersonProtect](#)

PRODUCT DESCRIPTION

Disposable KN95 Face Mask

Material: 4ply(non-woven+Hot air cotton+meltblown cloth+non-woven)

Color: White

Standard: China GB2626-2006

Application:Personal protection,PM2.5 protection

Filtering effect:95%

Product Features

- 1.Multi layer filter structure protection
- 2.Ear belt / head-mounted
- 3.Comfortable and breathable
- 4.Good protective effect
- 5.Comfortable soft elastic band
- 6.Soft sponge nose pad
- 7.Aluminum adjustable nose clip

General information - Contact us for more details if needed:

Supplied from multiple factories

All masks quality checked and CE certificates verified

CE certificate / Test Certificates / Photos can vary by manufacturer.

Packaging: Variable. Op-Bags 1/4 masks / Larger bags 10/20/50 masks.

Typical QTY in carton: 1000 masks.

Special packing configuration available upon request.

[Bioide-Type95_TechFile220520](#)

